


# Obsidian®

## MILLING BLOCKS



### THE LITHIUM SILICATE WITH EASY DELIVERY AND GREAT RESULTS

- Recommended for crowns, veneers, inlays and onlays
- Mills in under 20 minutes\*
- Complete library of shades that blend well with surrounding dentition
- Excellent long-term clinical performance



**SUPERIOR CHIP  
RESISTANCE**



**PATIENT-PLEASING  
RESULTS**



**WEAR-  
FRIENDLY**

Before



After

This patient presented with failing PFMs and recurrent decay on teeth #7 and #8. Obsidian® crowns were delivered in just one visit. To improve overall symmetry, same-day Obsidian veneers were placed on teeth #9 and #10.

“

*I've used CAD/CAM in my practice for the past 10 years, mostly to mill IPS e.max® crowns for same-day delivery. Since I started milling Obsidian blocks instead, I've been blown away by how esthetic in-office crowns can really be. The results are constantly impressive, rivaling those of any dental lab.*

”

– Christopher A. Petrush, DDS | Pleasant Hill, CA

For more information  
[glidewell.com/milling-blocks](http://glidewell.com/milling-blocks) ■ 844-949-7184



## BENEFITS

- Provide patients with lab-quality Obsidian® crowns in just one appointment
- Select from bonding or cementing delivery methods to suit each case
- Choose from a wide range of lifelike shades that blend well with existing dentition



Sold in 5-block packs

## SPECIFICATIONS

### MATERIAL

- Lithium silicate ceramic
- Flexural strength: 385 MPa

### IDEAL FOR


- Full-coverage crowns
- Inlays and onlays
- Veneers

## SHADES

- 12 VITA® Classical Shades: A1, A2, A3, A3.5, B1, B2, B3, C1, C2, C3, D2, D3
- 2 Ivoclar™ Bleach Shades: BL1, BL4

## EASY PREPARATION

### Anterior Full-Coverage Crown



\*Milling times vary based on prescription and mill used. Average milling time is shown.  
VITA is a registered trademark of VITA Zahnfabrik. Ivoclar is a registered trademark of Ivoclar Vivadent.

MKT-10576\_4 GD-313515-010820