

Foil trial package for ERKOFORM-RVE

Foil material	Ord.no. (regular contents)	Quantity	Thickn.	Applications (please see instructions in the thermoforming technique booklet and material card)
ERKOPLAST-O opaque-white	55 13 30 (45) 55 13 40 (40)	1 pc. 1 pc.	3.0 mm 4.0 mm	individual impression trays, upper jaw individual impression trays, lower jaw
ERKORIT clear	56 11 25 (50) 56 11 35 (40)	1 pc. 1 pc.	2.5 mm 3.5 mm	clear individual impression trays, upper jaw clear individual impression trays, lower jaw
ERKOLEN opaque-natural	51 11 10 (50)	2 pcs.	1.0 mm	temporary appliances, spacer for gel-type therapeutics, gingival dressings
ERKODUR clear	52 11 08 (60) 52 11 10 (50) 52 11 20 (30) 52 11 30 (20) 52 11 40 (15) 52 11 50 (12)	1 pc. 3 pcs. 3 pcs. 1 pc. 2 pcs. 2 pcs.	0.8 mm 1.0 mm 2.0 mm 3.0 mm 4.0 mm 5.0 mm	stabilization splints, retainers (Essix) stabilization splints, dressing plates occlusal splints, adjustment by reduction or addition (RESILIT-S) occlusal splints, adjustment by reduction or addition (RESILIT-S) occlusal splints, adjustment by reduction (OCCLUFORM) occlusal splints, adjustment by reduction (OCCLUFORM)
ERKOFLEX transparent	58 11 20 (12) 58 11 40 (6)	2 pc. 1 pc.	2.0 mm 4.0 mm	bracket transfer splints, base for PLAYSAFE sports mouthguards positioners, sports mouthguards, duplication
ERKOFLEX color	- (6)	1 pc.	4.0 mm	positioners, sports mouthguards
ERKOFLEX-freestyle confetti	58 11 85 (6)	1 pc.	4.0 mm	positioners, sports mouthguards
ERKOFLEX-bleach transparent	58 13 10 (20)	2 pcs.	1.0 mm	bleaching and fluoride trays
ERKOFLEX-95 opaque	58 95 40 (6)	1 pc.	4.0 mm	occlusal splints, single-layered sports mouthguards
ERKOLOC-pro hard/soft transparent	59 51 30 (10) 59 51 40 (10)	1 pc. 1 pc.	3.0 mm 4.0 mm	comfortable occlusal splints, adjustment by reduction or addition (RESILIT-S) comfortable occlusal splints, adjustment by reduction (OCCLUFORM)
Universal-UZF white	53 11 09 (50)	3 pcs.	0.1 mm	with ERKOFLEX (transparent) for better transparency

Sample/Material	Ord.no. (regular contents)	Quantity	Thickn.	Applications/Use
Cover templates	110 900 (8)	4 pcs.	-	to cover the filling granules when thermoforming soft materials such as ERKOFLEX, ERKOFLEX-95 and ERKOLOC/-pro
ERKOPOR	110 845 (40)	2 pcs.	3.0 mm	as a spacer when fabricating individual impression trays (app. 0.7 mm space per ERKOPOR-foamdisc after thermoforming)
ERKOGUM	110 844 (400 g)	app. 10 g	-	reusable blocking out material, silicone based. Especially suited for the protection of residual teeth, for blocking out undercuts and for placing artificial teeth for the fabrication of temporary appliances.
ISOLAC	624 050 (500 ml)	app. 50 ml	-	colorfree insulation, plaster against acrylics. For the fabrication of transparent splints out of ERKOLOC and ERKOFLEX, insulate plaster model prior to thermoforming.
ERKOSKIN	625 050 (50 ml)	app. 2,5 ml	-	as a spacer for bleaching trays and to relieve tension of the gingival margin for occlusal splints
LISKO sample	- (10)	1 pc. each assorted	-	LISKO-S unique for the finishing of soft materials and to finish the edges of all thermoforming materials. LISKO white (fine grain), LISKO grey (middle grain) and LISKO brown (coarse grain), as above, to finish soft relinings and to matt metal.
Tray handle out of plastic	222 901 (100)	1 pc.	-	to reinforce lower individual trays
Tray handle out of aluminium, short	222 857 (50)	1 pc.	-	to melt into individual trays
Tray handle out of aluminium, long	222 858 (50)	1 pc.	-	to melt into individual trays (larger handle)
Blocking out wax	725 080 (50 g)	app. 4 g	-	colorfree, high fusing wax to block out bubbles in the plaster
Magnet	250 151 (6)	1 pc.	-	to gather filling granules (high-grade steel, magnetic)
Brush red	188 010 (3)	1 pc.	-	to place and condense filling granules around the model